

# Certificate of Analysis

Dimethyl Sulfoxide, Cell culture tested, >99.7% , 0.2um filtered

Cat No : D2680

Lot No 05022208

TEST	Specification	Result
Appearance	Fine White to Off White Powder	PASS
Identification	Confirmed	PASS
Purity (Dried based by GC)	$\geq 99.0 \%$	$\geq 99.9 \%$
Density at 25°C	$\geq 1.095$	1.0953
Titration acid	$\leq 0.0005 \text{ meq/g}$	0.0001 meq/g
Residue after Evaporation	$\leq 0.0005 \%$	< 0.001%
Water	$\leq 0.1 \text{ (KF)}$	< 0.05%

This document has been produced electronically and is valid without signature

Date of Print : 5. 4. 2022